

Numer Katalogowy	Nazwa Produktu	Wielkość Opakowania
M0491S	Q5 High-Fidelity DNA Polymerase	100 units
M0491L		500 units
M0493S	Q5 Hot Start High-Fidelity DNA Polymerase	100 units
M0493L		500 units
M0492S	Q5 High-Fidelity 2X Master Mix	100 rxns (50 µl vol)
M0492L		500 rxns (50 µl vol)
M0494S	Q5 Hot Start High-Fidelity 2X Master Mix	100 rxns (2x1.25ml)
M0494L		500 rxns (10x1.25ml)
M0494X		500 rxns (1x12.5ml)
E0555L	Q5 High-Fidelity PCR Kit	200 rxns
E0555S		50 rxns
M0515S	Q5U Hot Start High-Fidelity DNA Polymerase	100 units
M0515L		500 units
M0500S	Q5 Blood Direct 2X Master Mix	100 rxns (1x1ml)
M0500L		500 rxns (5x1ml)
M0578L	Q5 Quick-Load High-Fidelity 2X Master Mix	500 rxns (10x1.25ml)
M0580L	Q5 Quick-Load Hot Start High-Fidelity 2X Master Mix	500 rxns (10x1.25ml)
M2499S	Q5-XT Hot Start High-Fidelity 2X Master Mix	100 rxns (2x1.25ml)
M2499L		500 rxns (10x1.25ml)
M0544S	NEBNext Ultra II Q5 Master Mix	50 rxns
M0544L		250 rxns
M0544X		500 rxns
M0597S	NEBNext Q5U Master Mix	50 rxns
M0597L		250 rxns